

**SEE** it  
EVEN **BETTER**™

**LED**  
DayLite™



Adding Light INCREASES depth of field,  
CONTRAST AND OVERALL visibility.

SEE IT EVEN BETTER WITH

**DESIGNS  
FOR VISION**®

**LED**  
DayLite™

**UNTETHERED, LIGHTWEIGHT AND PORTABLE**

### **NEW TWIN BEAM™ HEADSET**

Optically collimated for even brighter light, the advanced LED DayLite™ Twin Beam™ Headset provides up to 13,000 foot-candles (139,930 lux) of WHITE light (5800°K). This lightweight and untethered light system is comparable in light output to the heavier, more costly xenon. Moreover, LEDs do not degrade and can be expected to provide illumination for thousands of hours!

Simply look for the Blue Cable™ and

**JUST SEE it!**™

with our 45 Day, NO Obligation Trial.  
Simply contact us for an appointment.

Call Toll Free **1.800.345.4009**

Or Visit Us @ [www.designsforvision.com](http://www.designsforvision.com)

LED DayLite™ Twin Beam™ Headset	
Light Intensity	13,000 foot-candles (139,930 lux) at 16"
Spot Size	3 inches at 16" working distance
Correlated Color Temperature	5800° K
Light Output	User Adjustable, 0-100%
Battery Pack w/belt clip Run Time	Lithium Ion (2 supplied) 4.5 hours per pack, 9 hours with Power Link Cable
Recharge Time	3 Hours
Headset Weight	7.0 oz.
Headset Warranty	Lifetime



**DESIGNS  
FOR VISION®**

**LED**  
*DayLite™*

**UNTETHERED, LIGHTWEIGHT AND PORTABLE,** the advanced LED DayLite™ provides bright WHITE Light for all of your operatory needs. The Blue Cable™ connects to a battery pack on your belt, is shielded for your protection and contains a Tellurium Copper Alloy Heat Sink to provide a cooler, longer lasting LED. Available in your choice of a Small Optically Focused LED DayLite™ that clips onto your telescopes (or is attached to a Lightweight Headband) or the NEW TwinBeam™ LED DayLite™ Headset, optically collimated to provide light bright enough to rival any xenon.



**JUST SEE it!™**

with our 45 Day, NO Obligation Trial.  
Simply contact us for an appointment.

Call Toll Free **1.800.345.4009**

Or Visit Us @ [www.designsforvision.com](http://www.designsforvision.com)

**DESIGNS FOR VISION™**

760 Koehler Avenue | Ronkonkoma, NY 11779 | 1.800.345.4009 | 1.631.585.3300 | [www.designsforvision.com](http://www.designsforvision.com)